




AED Readiness

*Peace of mind with
Arch™ AED Program
Management and
Medical Direction*



An AED Superstore Product

Easily Manage Your AED Program



AED Regulations vary by state and may include any of the following requirements:

- Ongoing physician oversight
- Physician approved policies and procedures manual
- CPR training and certification
- Registering your AED with the local dispatch center or EMS agency
- Performing and documenting readiness checks every 30 days
- Submission of reports after the AED is used

While not complex, these requirements can be time consuming if managed manually. Arch is your turn-key solution.



Arch provides comprehensive AED management

- Monthly reminders to check AEDs with late alerts.
- Notification when electrode pads and batteries are about to expire and access to our online store for easy reordering.
- Online access to your personalized Policy and Procedures Manual.
- Automatic registration with your local 911 call center & local EMS (if required).
- Tracking of CPR certification expiration dates for all your trained responders.
- Physician oversight specific to your state's requirements.
- Easy access to a support team to answer questions or discuss concerns.
- Online support ticketing to report issues with your AED.
- Online AED usage reporting and support:
 - Next day shipping of same make and model loaner AED.
 - Replacement AED electrode pads.
 - All shipping labels.
 - Data download from your AED.
 - Physician review of data and post AED event report.
 - Submission of required regulatory reports as needed per your location's requirements.

Global Administrator Homepage with Dashboard

ORGANIZATION DETAILS:

AED UNITS LOCATIONS RESPONDERS & USERS DASHBOARD

AED Equipment Status	
All AEDs	48
Expired AEDs	6
Expiring AEDs	5
Up-to-date AEDs	37

AED Readiness Status	
All AEDs	48
Overdue Readiness	17
Expiring Readiness	22
Up-to-date Readiness	9

Trained Responders	
All Trained Responders	27
Overdue	9
Expiring within 30 days	1
Up-to-date	13
Unknown	4

Homepage Tab Views

ORGANIZATION DETAILS:

AED UNITS LOCATIONS RESPONDERS & USERS DASHBOARD

Click the icon under the Readiness column to perform a readiness check. Click the icon under the Expiration column to update expiration.

Model	Location	Serial Number	Placement	Readiness	Expiration
Cardiac Science Powerheart G3	Corporate Headquarters	999999	to the left of the bathroom	🔴	🔴
main battery expired					
Cardiac Science Powerheart G3	Corporate Headquarters	999999	Lunch Room	🟡	🟢
Cardiac Science Powerheart G3 Pro	Corporate Headquarters	999997	Guard Shack	🔴	🟢
Defitech LibLife	Corporate Headquarters	1010000	By Bathrooms	🟢	🟢
Defitech LibLife VIEW	Corporate Headquarters	4000000	Health Center	🔴	🟢
HeartSine Samaritan PAD 350P	Corporate Headquarters	100000000	Security Desk	🔴	🔴

ORGANIZATION DETAILS:

AED UNITS LOCATIONS RESPONDERS & USERS DASHBOARD

Name	Overall Loc Status	#AED Units	Readiness Alerts	Exp Alerts	Exp Warnings
Corporate Headquarters	🔴	15	13	2	0
Demo - Blismard	🔴	5	1	1	0
Demo - Boston	🟡	3	0	0	2
Demo - Little Rock	🟢	0	0	0	0
Demo - Los Angeles	🟢	4	0	0	0
Demo - Madison	🟢	4	0	0	0

AED Units Tab

ORGANIZATION DETAILS:

AED UNITS LOCATIONS RESPONDERS & USERS DASHBOARD

Name	Location	Email	Roles	Conv Status
John Smith			Responder	🟢
Kayla Chan		kchan@demoaccount.com	Responder	🔴
Cecilia Arnoldson	Demo - Blismard	carnoldson@demoaccount.com	Local Admin, Responder	🔴
Michael Peters	Demo - Blismard	mpeters@demoaccount.com	Responder	🔴
Beth Marks	Demo - Blismard	bmarks@demoaccount.com	Responder	🟡
John Jacobs	Demo - Blismard	jjacobs@demoaccount.com	Responder	🟡

Responders & Users Tab

Locations Tab

ORGANIZATION DETAILS:

AED UNITS LOCATIONS RESPONDERS & USERS DASHBOARD

AED Equipment Status	
All AEDs	48
Expired AEDs	6
Expiring AEDs	5
Up-to-date AEDs	37


AED Readiness Status	
All AEDs	48
Overdue Readiness	17
Expiring Readiness	22
Up-to-date Readiness	9

Trained Responders	
All Trained Responders	27
Overdue	9
Expiring within 30 days	1
Up-to-date	14
Unknown	4

Dashboard Tab

How to Record a Readiness Check

AED Readiness through the Device Page



Alternate way to enter check is at the device level. From home page with map, click on AED serial number under Organization Details

Access to view QR Code for an easy way to record AED readiness check

*Spare battery, spare set of electrodes, or pediatric electrodes expiration dates will also be displayed here. Manage the list by selecting "Set New Expiration Dates"



QR Codes

Quick-Response Codes are an easy way to check the AED with your smartphone

